POSTER PRESENTATION





Evidence for higher sensitivity of recombinant Tri a 36 compared to omega-5-gliadin for diagnosis of wheat food allergy

A Baar^{1*}, S Pahr¹, C Constantin², S Giavi³, M Alkistib³, NG Papadopoulos³, C Ebner⁴, A Mari⁵, S Vrtala¹, R Valenta⁶

From Food Allergy and Anaphylaxis Meeting (FAAM 2013) Nice, France. 7-9 February 2013

Background

Wheat is one of the most important food allergen sources. Using natural wheat allergen extracts for serological diagnosis of wheat-induced food allergy false positive test results are frequently obtained, in particular in grass pollen allergic patients. Therefore, Tri a 19, an omega-5-gliadin, which is known as a major allergen in wheat dependent exercise induced anaphylaxis (WDEIA) and in wheat food allergy in children, is widely used for the serological diagnosis of wheat-induced food allergy.

Methods

We have recently characterized Tri a 36, a low molecular weight glutenin as new major wheat food allergen. Here we compared recombinant Tri a 36 with ImmunoCAPs containing omega-5-gliadin for the detection of specific IgE antibodies in a population of wheat food allergic children (n=23) and grass pollen allergic patients (n=21). Wheat food allergic children had a clear clinical history of wheat food allergy with symptoms clearly attributable to wheat ingestion. Grass pollen allergic patients suffered from grass pollen-induced respiratory allergy but regularly ate wheat products without any clinical symptoms. Sera from both populations were tested by ImmunoCAPs containing wheat extract or omega-5-gliadin and by IgE ELISA to rTri a 36.

Results

Using wheat extract-based ImmunoCAPs all but one of the wheat food allergic patients showed allergen-specific IgE levels >0.35 kUA/L but 17 out of the 21 grass pollen allergic patients gave false positive test results. Eleven (i. e.,48%) of the wheat food allergic patients showed >0.35 kUA/L IgE against omega-5-gliadin and two of the grass pollen allergic patients gave false positive test results. Using rTri a 36, 14 (i.e., 61%) of the wheat food allergic patients were diagnosed and only one grass pollen allergic patient gave a false positive test result. Each of the sera from the wheat food allergic patients with IgE reactivity to omega-5-gliadin also reacted with rTri a 36.

Conclusion

Our results indicate that rTri a 36 has higher sensitivity and specificity than omega-5-gliadin for the diagnosis of wheat food allergy.

Disclosure of interest

None declared.

Acknowledgement

This study was supported by the Christian Doppler Research Association, Vienna Austria, by a research grant from Phadia/Thermofisher, Uppsala, Sweden and by European Commission's Seventh Framework Programme under grant agreement No. 261357.

Full list of author information is available at the end of the article



© 2013 Baar et al; licensee BioMed Central Ltd. This is an Open Access article distributed under the terms of the Creative Commons Attribution License (http://creativecommons.org/licenses/by/2.0), which permits unrestricted use, distribution, and reproduction in any medium, provided the original work is properly cited.

¹Department of Pathophysiology and Allergy Research, Christian Doppler Laboratory for the Development of Allergen Chips, Medical University of Vienna, Austria

Author details

¹Department of Pathophysiology and Allergy Research, Christian Doppler Laboratory for the Development of Allergen Chips, Medical University of Vienna, Austria. ²Department of Pathophysiology and Allergy Research, Medical University Vienna, Vienna, Austria. ³Allergy and Immunology Research Centre, University of Athens, Athens, Greece. ⁴Ambulatory for Allergy and Clinical Immunology, Ambulatory for Allergy and Clinical Immunology, Vienna, Austria. ⁵Center for Molecular Allergology, IDI-IRCCS, Rome, Italy. ⁶Department of Pathophysiology and Allergy Research, Christian Doppler Laboratory for Allergy Research, Medical University Vienna, Vienna, Austria.

Published: 25 July 2013

doi:10.1186/2045-7022-3-S3-P171

Cite this article as: Baar *et al.*: **Evidence for higher sensitivity of** recombinant **Tri a 36 compared to omega-5-gliadin for diagnosis of wheat food allergy.** *Clinical and Translational Allergy* 2013 **3**(Suppl 3):P171.

Submit your next manuscript to BioMed Central and take full advantage of:

- Convenient online submission
- Thorough peer review
- No space constraints or color figure charges
- Immediate publication on acceptance
- Inclusion in PubMed, CAS, Scopus and Google Scholar
- Research which is freely available for redistribution

Submit your manuscript at www.biomedcentral.com/submit

BioMed Central